

Product Datasheet

CABYR peptide (orb218475)

Description	The peptide is used to block anti-CBP86 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CABYR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215119.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D424, O75952
Entrez	71132, 26256
Expiration Date	12 months from date of receipt.