

Product Datasheet

RANBP6 peptide (orb218481)

Description	The peptide is used to block anti-RANBP6 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RANBP6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215125.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BIV3, O60518
Entrez	26953, 240614
Expiration Date	12 months from date of receipt.