

Product Datasheet

ATP5S peptide (orb218486)

Description	The peptide is used to block anti-ATP5S antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ATP5S
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215130.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99766, Q9CRA7, Q5XIM4
Entrez	27109, 68055, 362749
Expiration Date	12 months from date of receipt.