

Product Datasheet

IL37 peptide (orb218488)

Description	IL37 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	IL37
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215132.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NZH6
Entrez	27178
Expiration Date	12 months from date of receipt.