

Product Datasheet

KCNMB4 peptide (orb218496)

Description	The peptide is used to block anti-KCNMB4 antibody reactivity.
Reactivity	Human, Porcine, Sheep
Tested Applications	Blocking
Target	KCNMB4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215140.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q86W47
Entrez	27345
Expiration Date	12 months from date of receipt.