

Product Datasheet

HIPK2 peptide (orb218499)

Description	The peptide is used to block anti-HIPK2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	HIPK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215143.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9QZR5, Q9H2X6
Entrez	28996, 15258
Expiration Date	12 months from date of receipt.