

Product Datasheet

SENP1 peptide (orb218506)

Description	The peptide is used to block anti-SENP1 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	SENP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215150.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9P0U3
Entrez	29843
Expiration Date	12 months from date of receipt.