

Product Datasheet

GPR171 peptide (orb218507)

Description	The peptide is used to block anti-GPR171 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GPR171
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215151.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BG55 , O14626
Entrez	29909 , 229323
Expiration Date	12 months from date of receipt.