

Product Datasheet

SLC6A18 peptide (orb218515)

Description	The peptide is used to block anti-XTRP2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC6A18
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215159.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q96N87 , O88576 , Q62687
Entrez	22598 , 29323 , 348932
Expiration Date	12 months from date of receipt.