

Product Datasheet

NOX4 peptide (orb218516)

Description	NOX4 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NOX4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215160.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q924V1 , Q9NPH5 , Q9JHI8
Entrez	50490 , 85431 , 50507
Expiration Date	12 months from date of receipt.