

Product Datasheet

ANGPTL4 peptide (orb218522)

Description	ANGPTL4 peptide
Reactivity	Human
Tested Applications	Blocking
Target	ANGPTL4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215166.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9BY76
Entrez	51129
Expiration Date	12 months from date of receipt.