

Product Datasheet

A4GNT peptide (orb218523)

Description	The peptide is used to block anti-A4GNT antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	A4GNT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215167.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UNA3
Entrez	51146
Expiration Date	12 months from date of receipt.