

Product Datasheet

SLC22A17 peptide (orb218526)

Description	The peptide is used to block anti-SLC22A17 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC22A17
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215170.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D9E0 , Q9P290 , Q8WUG5
Entrez	51310 , 59049 , 305886
Expiration Date	12 months from date of receipt.