

Product Datasheet

UBR5 peptide (orb218529)

Description	The peptide is used to block anti-UBR5 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Zebrafish
Tested Applications	Blocking
Target	UBR5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215173.
Purity	> 85%
Source	Synthetic
Uniprot ID	O95071, Q80TP3
Entrez	70790, 51366
Expiration Date	12 months from date of receipt.