

Product Datasheet

ANGPT4 peptide (orb218530)

Description	The peptide is used to block anti-Angiopoietin-4 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ANGPT4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215174.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Y264, Q9WVH6
Entrez	51378, 11602
Expiration Date	12 months from date of receipt.