

Product Datasheet

WNT16 peptide (orb218531)

Description	The peptide is used to block anti-WNT16 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	WNT16
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215175.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UBV4
Entrez	51384
Expiration Date	12 months from date of receipt.