

Product Datasheet

IP6K2 peptide (orb218534)

Description	The peptide is used to block anti-IP6K2 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	IP6K2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215178.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UHH9
Entrez	51447
Expiration Date	12 months from date of receipt.