

Product Datasheet

HDAC7 peptide (orb218537)

Description	The peptide is used to block anti-Histone Deacetylase 7 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	HDAC7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215181.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8WUI4 , Q99P96 , Q8C2B3
Entrez	56233 , 51564
Expiration Date	12 months from date of receipt.