

Product Datasheet

DNAJB11 peptide (orb218546)

Description	The peptide is used to block anti-DNAJB11 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DNAJB11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215190.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6TUG0, Q99KV1, Q9UBS4
Entrez	51726, 360734, 67838
Expiration Date	12 months from date of receipt.