

Product Datasheet

ARID4B peptide (orb218547)

Description	The peptide is used to block anti-ARID4B antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	ARID4B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215191.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JKB5, Q4LE39
Entrez	84481, 51742
Expiration Date	12 months from date of receipt.