

Product Datasheet

SLC6A15 peptide (orb218564)

Description	The peptide is used to block anti-SBAT1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SLC6A15
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215208.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H2J7, Q08469, Q8BG16
Entrez	103098, 282712, 55117
Expiration Date	12 months from date of receipt.