

Product Datasheet

SLC39A9 peptide (orb218570)

Description	The peptide is used to block anti-ZIP9 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC39A9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215214.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BFU1, Q9NUM3, Q3KR82
Entrez	328133, 314275, 55334
Expiration Date	12 months from date of receipt.