

Product Datasheet

CXCR7 peptide (orb218591)

Description	The peptide is used to block anti-GPR159 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CXCR7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215235.
Purity	> 85%
Source	Synthetic
Uniprot ID	P56485, P25106, O89039
Entrez	57007, 84348, 12778
Expiration Date	12 months from date of receipt.