

Product Datasheet

CCNL1 peptide (orb218593)

Description	The peptide is used to block anti-Cyclin L1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	CCNL1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215237.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q52KE7, Q9UK58, Q9R1Q2
Entrez	114121, 57018, 100909712, 56706
Expiration Date	12 months from date of receipt.