

Product Datasheet

CHRNA10 peptide (orb218594)

Description	The peptide is used to block anti-CHRNA10 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	CHRNA10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215238.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9GZZ6
Entrez	57053
Expiration Date	12 months from date of receipt.