

Product Datasheet

SENP7 peptide (orb218598)

Description	The peptide is used to block anti-SENP7 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	SENP7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215242.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9BQF6, Q8BUH8, D3ZF42
Entrez	57337, 66315, 288167
Expiration Date	12 months from date of receipt.