

Product Datasheet

COL20A1 peptide (orb218600)

Description	The peptide is used to block anti-Collagen 20 alpha 1 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	COL20A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215244.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9P218
Entrez	57642
Expiration Date	12 months from date of receipt.