

Product Datasheet

PKNOX2 peptide (orb218610)

Description	The peptide is used to block anti-PKNOX2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PKNOX2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215254.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BG99 , Q96KN3
Entrez	208076 , 63876
Expiration Date	12 months from date of receipt.