

Product Datasheet

ANAPC1 peptide (orb218622)

Description	The peptide is used to block anti-ANAPC1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ANAPC1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215266.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H1A4, P53995
Entrez	17222, 64682
Expiration Date	12 months from date of receipt.