

Product Datasheet

IPPK peptide (orb218624)

Description	The peptide is used to block anti-IPPK antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	IPPK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215268.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5PXE9, Q9H8X2, Q6P1C1
Entrez	75678, 306808, 64768
Expiration Date	12 months from date of receipt.