

Product Datasheet

CDCP1 peptide (orb218628)

Description	The peptide is used to block anti-CD318 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CDCP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215272.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5U462, Q9H5V8
Entrez	109332, 64866
Expiration Date	12 months from date of receipt.