

Product Datasheet

FTO peptide (orb218639)

Description	The peptide is used to block anti-FTO antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	FTO
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215283.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BGW1, Q2A121, Q9C0B1
Entrez	291905, 26383, 79068
Expiration Date	12 months from date of receipt.