

Product Datasheet

COASY peptide (orb218658)

Description	The peptide is used to block anti-COASY antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	COASY
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215302.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9DBL7, Q13057
Entrez	71743, 80347
Expiration Date	12 months from date of receipt.