

Product Datasheet

SLC25A31 peptide (orb218665)

Description	The peptide is used to block anti-ANT4 antibody reactivity.
Reactivity	Human, Primate
Tested Applications	Blocking
Target	SLC25A31
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215309.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H0C2
Entrez	83447
Expiration Date	12 months from date of receipt.