

Product Datasheet

SLC4A9 peptide (orb218668)

Description	SLC4A9 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC4A9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215312.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q96Q91 , Q8K4V2
Entrez	83697 , 266612
Expiration Date	12 months from date of receipt.