

Product Datasheet

SYT16 peptide (orb218669)

Description	The peptide is used to block anti-SYT16 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	SYT16
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215313.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q17RD7
Entrez	83851
Expiration Date	12 months from date of receipt.