

Product Datasheet

HORMAD1 peptide (orb218674)

Description	The peptide is used to block anti-HORMAD1 antibody reactivity.
Reactivity	Human, Mouse, Primate
Tested Applications	Blocking
Target	HORMAD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215318.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q86X24
Entrez	84072
Expiration Date	12 months from date of receipt.