

Product Datasheet

RAB11FIP4 peptide (orb218677)

Description	The peptide is used to block anti-RAB11FIP4 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RAB11FIP4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215321.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q86YS3, Q8BQP8
Entrez	84440, 268451
Expiration Date	12 months from date of receipt.