

Product Datasheet

GPR174 peptide (orb218680)

Description	The peptide is used to block anti-GPR174 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GPR174
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215324.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q3U507 , Q9BXC1
Entrez	213439 , 84636
Expiration Date	12 months from date of receipt.