

Product Datasheet

GNG8 peptide (orb218700)

Description	The peptide is used to block anti-GNG8 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep
Tested Applications	Blocking
Target	GNG8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215344.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UK08 , P63078 , P63077
Entrez	14709 , 245986 , 94235
Expiration Date	12 months from date of receipt.