

Product Datasheet

HIST1H4A peptide (orb218711)

Description	The peptide is used to block anti-Histone H4
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	HIST1H4A; HIST1H4B; HIST1H4C; HIST1H4D; HIST1H4E; HIST1H4F; HIST1H4H; HIST1H4I; HIST1H4J; HIST1H4K; HIST1H4L; HIST2H4A; HIST2H4B; HIST4H4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215355.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62806, P62804, P62805
Entrez	319161, 320332, 100041230, 326620, 8368, 102557184, 8361, 680097, 554313, 8294, 8365, 500351, 102548682, 97122, 8364, 100912290, 502913, 8359, 8370, 100912418, 102641229, 291152, 319160, 8366, 100360950, 295277, 326619, 64627, 121504, 319158, 319159, 8363, 69386, 102551184, 8362, 100912564, 319155, 8367, 319157, 8360, 319156
Expiration Date	12 months from date of receipt.