

Product Datasheet

SENp8 peptide (orb218713)

Description	The peptide is used to block anti-SENp8 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SENp8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215357.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5FVJ8, Q96LD8, Q9D2Z4
Entrez	123228, 71599, 315723
Expiration Date	12 months from date of receipt.