

Product Datasheet

Apolipoprotein AIBP peptide (orb218719)

Description	Apolipoprotein AIBP peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	APOA1BP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215363.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8K4Z3 , B0BNM1 , Q8NCW5
Entrez	295229 , 128240 , 246703
Expiration Date	12 months from date of receipt.