

Product Datasheet

MYO3B peptide (orb218729)

Description	The peptide is used to block anti-MYO3B antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MYO3B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215373.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q1EG27 , Q8WXR4
Entrez	140469 , 329421
Expiration Date	12 months from date of receipt.