

Product Datasheet

CCBE1 peptide (orb218737)

Description	CCBE1 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CCBE1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215381.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q3MI99, Q6UXH8
Entrez	320924, 147372
Expiration Date	12 months from date of receipt.