

Product Datasheet

CBLN2 peptide (orb218738)

Description	The peptide is used to block anti-Cerebellin 2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	CBLN2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215382.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8IUUK8, P98087, Q8BGU2
Entrez	147381, 12405, 291388
Expiration Date	12 months from date of receipt.