

Product Datasheet

DSG4 peptide (orb218739)

Description	The peptide is used to block anti-Desmoglein 4 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DSG4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215383.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6W3B0, Q7TMD7, Q86SJ6
Entrez	147409, 291754, 16769
Expiration Date	12 months from date of receipt.