

Product Datasheet

TCEB3C peptide (orb218747)

Description	The peptide is used to block anti-TCEB3C antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TCEB3C
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215391.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8NG57
Entrez	102725042, 101930171, 162699, 101930165
Expiration Date	12 months from date of receipt.