

Product Datasheet

SLC16A13 peptide (orb218754)

Description	The peptide is used to block anti-MCT13 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC16A13
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215398.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7RTY0 , Q66HE2 , Q8CE94
Entrez	287451 , 201232 , 69309
Expiration Date	12 months from date of receipt.