

Product Datasheet

SEN5 peptide (orb218756)

Description	The peptide is used to block anti-SEN5 antibody reactivity.
Reactivity	Human, Porcine, Primate, Rat
Tested Applications	Blocking
Target	SEN5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215400.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q96HI0
Entrez	205564
Expiration Date	12 months from date of receipt.